


Salmon water to be collected in 1.25 L (1 qt) bottles, to discharge the scale, always labeled in μ g/liter.

The three studies that are included in our meta-analysis have a 90% confidence interval of 0.0001, and the three studies that are excluded have a 95% confidence interval of 0.0001. We did not find any publication bias. No significant difference was found between the 90% confidence interval of 0.0001 and the 95% confidence interval of 0.0001. The three studies that are included in our meta-analysis have a 90% confidence interval of 0.0001, and the three studies that are excluded have a 95% confidence interval of 0.0001. We did not find any publication bias. No significant difference was found between the 90% confidence interval of 0.0001 and the 95% confidence interval of 0.0001.



 **LABC** partner
Office 3K4,
Redwither Business Tower,
Wrexham, LL13 9XT
Tel: 01978 356500
www.bjreprintarchitectural.com
© copyright all rights reserved